

### First Aid



### The +

Ease of use  
Safety: medical device

#### PURCHASE PACKAGING

| Ref.  | Size   | Inner | Carton |
|-------|--------|-------|--------|
| 60115 | 500 ml |       | 12     |

#### SALES PACKAGING



Multiple

#### DESCRIPTION

Ophthalmic solution - accidents involving foreign bodies  
Sterile saline solution of sodium chloride (0.9%) in case of solid particles (dust, chips, metal splinters...) projections or liquids (chemicals such as oils or solvents).

500 ml bottle with an ergonomic eyepiece (not supplied with a cap anymore).

200 ml bottle - to be used individually or as a refill on eye wash stations and eye wash boxes.

Rinsing time: approx. 2 min.

#### SECTORS

Mining

Heavy & process industries

Mechanical constructions

Chemical & pharmaceutical industries

Light Industries

Regional/public authorities

Infrastructure, building industries

## TECHNICAL FEATURES

|            |                   |
|------------|-------------------|
| Shape      | First Aid         |
| Dimensions | 60 x 230 x 290 mm |
| Color      | White             |
| Color 2    | Green             |

## STANDARD(S)

This equipment conforms to the personal protective equipment model covered by the EC type-examination certificate **DGM-272**

Delivered by Presafe (0543) Tuborg Parkvej 8 DK-2900 Hellerup Denmark



DIRECTIVE 93/42/EEC ON  
MEDICAL