

Ref. MO60111

First Aid



The +

Ease of use
Safety: Medical device

PURCHASE PACKAGING

| Ref. | Size | Inner | Carton |
|-------|--------|-------|--------|
| 60111 | 200 ml | | 10 |

SALES PACKAGING

>1

Multiple

DESCRIPTION

PLUM bottle 200 ml pH Neutral - neutralizing solution
Sterile phosphate buffer solution (4.9%) neutralizing acids and bases. Very quickly brings the pH of the eye back to a normal value of 7.4.
200 ml bottle with an ergonomic eyepiece holder (not supplied with a cap anymore).
200 ml bottle to be used individually or as a refill on eyewash stations and eye wash sets.
Rinsing time: 2 min

SECTORS

- Mining
- Heavy & process industries
- Mechanicals constructions
- Chemical & pharmaceutical industries
- Light Industries
- Infrastructure, building industry

TECHNICAL FEATURES

| | | | |
|------------|-------------|---------|------|
| Shape | First Aid | Color 2 | Blue |
| Weight | 200 g | | |
| Dimensions | 165 x 55 mm | | |
| Color | White | | |

STANDARD(S)

This equipment conforms to the personal protective equipment model covered by the EC type-examination certificate **DGM-272**

Delivered by Presafe (0543) Tuborg Parkvej 8 DK-2900 Hellerup Denmark



DIRECTIVE 93/42/EEC ON
MEDICAL