

Ref. MO60116

First Aid



The +

Ease of use
Safety: medical device

PURCHASE PACKAGING

Référence	Taille	Sachet	Carton
60116			6

DESCRIPTION

500 ml pH Neutral DUO

Bottle 500 ml DUO Sterile neutralizing solution of phosphate buffer (4.9%) - accidents involving acidic and basic chemicals

Sterile phosphate buffer solution (4.9%) neutralizing acids and bases. Very quickly brings the pH of the eye back to a normal value of 7.4.

500 ml DUO bottle equipped with a double eye cup for simultaneous rinsing of both eyes.

500 ml DUO bottle to be used individually or as a refill for eye wash sets.

Rinsing time: approx. 2 min.

SECTORS

Heavy & process industries

Mechanicals constructions

Chemical & pharmaceutical industries

Light Industries

Regional/public authorities

Infrastructure, building industry

Construction - finishing work

TECHNICAL FEATURES

Color	White
Color 2	Blue
Dimensions	220 x 75 mm

STANDARD(S)

This equipment conforms to the personal protective equipment model covered by the EC type-examination certificate **DGM-272**

Delivered by Presafe (0543) Tuborg Parkvej 8 DK-2900 Hellerup Denmark



DIRECTIVE 93/42/EEC ON MEDICAL

RETAILER STAMP

MAJ 04/03/2021