

Ref. MO60110

First Aid



The +

Ease of use
Safety: medical device

PURCHASE PACKAGING

Ref.	Size	Inner	Carton
60110	200 ml		10

SALES PACKAGING

>1

Multiple

DESCRIPTION

Ophthalmic solution - accidents involving foreign bodies
Sterile saline solution of sodium chloride (0.9%) in case of solid particles (dust, chips, metal splinters...) projection or liquids (chemicals such as oils or solvents).

200 ml bottle with an ergonomic eyepiece a(not supplied with a cap anymore).

200 ml bottle to be used individually or as a refill on eye wash stations and eye wash boxes.

Rinsing time: approx. 2 min.

SECTORS

Heavy & process industries

Chemical & pharmaceutical industries

Light Industries

Infrastructure, building industry

TECHNICAL FEATURES

Shape	First Aid	Color 2	Green
Weight	200g		
Dimensions	165 x 55 mm		
Color	White		

STANDARD(S)

This equipment conforms to the personal protective equipment model covered by the EC type-examination certificate **DGM-272**

Delivered by Presafe (0543) Tuborg Parkvej 8 DK-2900 Hellerup Denmark



DIRECTIVE 93/42/EEC ON
MEDICAL