

## **TECHNICAL DATASHEET**

# **PLUM eyewash Station**

Ref. MO60120

First Aid





#### **PURCHASE PACKAGING**

RéférenceTailleSachetCarton60120500ml1

#### **DESCRIPTION**

#### PLUM eyewash Station

Open wall station with 2 Plum 500 ml eye wash bottles. Suitable for dust-free work areas where there is a risk of foreign bodies entering the eyes.

Sterile saline solution of sodium chloride (0.9%) in case of solid particles (dust, shavings, metal splinters...) projection or liquids (chemicals such as oils or solvents).

Dimensions: L 290 x W 228 x H 80 mm

Rinsing time: approx. 10 min

### **TECHNICAL FEATURES**

**Dimensions** 

60 x 230 x 290 mm

## STANDARD(S)

This equipment conforms to the personal protective equipment model covered by the EC type-examination certificate **DGM-272** Delivered by Presafe (0543) Tuborg Parkvej 8 DK-2900 Hellerup Denmark



**DIRECTIVE 93/42/EEC ON MEDICAL** 

RETAILER STAMP

MAJ 04/03/2021